

## AccuRx in virtual wards

### Clinician Guide

1. Patient will be registered onto the virtual ward as a temporary or walk in patient
2. You can send the patient accuRx Florey templates in the usual way by clicking on the message icon on the accuRx toolbar



3. Click on Florey and choose the Covid-19 monitoring questionnaire. You can choose what free-text is sent to accompany your Florey template. You can also set up the scheduler to send a daily text message by clicking on the calendar icon. Please note you can not send a new message until the last text has been replied to.

**POWER, Michelle (Mrs)**  
NHS No. 289-646-0683    15-Sep-1926 (94y)    Female

**Consent given for communication by SMS text messaging (04-01-18)**  
Declined consent for SMS text messaging (04-01-18)

Michelle's number  
07894561230

Template | Pathway | Florey | Video

Click to search questionnaires, e.g. Asthma

- Asthma ACT Questionnaire
- Asthma RCP Questionnaire
- COPD Questionnaire
- COPD Questionnaire - CAT
- COPD Questionnaire - MRC
- COVID-19 Remote Monitoring Questionnaire
- COVID-19 Triage Questionnaire
- Flu Vaccination Invite
- Lifestyle surveys (requires covering SMS costs)**
  - Alcohol Consumption Questionnaire - AUDIT-C
  - BMI Questionnaire
  - Smoking Questionnaire

Click here  
to change  
the date



Full information on how to set this up is [HERE<sup>1</sup>](#)

4. Patient replies will return to the accuRx inbox. If you are a member of the remote monitoring 'team' you will see a red notification (below)



5. Read and save patient replies to the patient record from within the accuRx inbox

More information on reading and actioning Florey replies is [HERE<sup>2</sup>](#)

6. The following information will be entered into the clinical record (EMIS or Systmone). The clinical readings will be recorded as SNOMED codes

@COVID-19 Team patient NAME has completed their COVID-19 Remote Monitoring Questionnaire.

Questionnaire: COVID-19 Remote Monitoring Questionnaire

DoB: 00-00-0000

Today vs yesterday: Better/Same/Worse

More breathless today: Yes/No

Full sentence without having to take a breath: Yes/No

Access to equipment: Thermometer, Blood pressure monitor, Fingertip oxygen monitor

Temperature (°C): 00.0

SYS BP (mmHg): 000

DIA BP (mmHg): 00

SpO2%: 00

Pulse (bpm): 00

7. If the patient is discharged from the service, or no longer needs the scheduled remote monitoring texts click [HERE<sup>3</sup>](#) to find out how to cancel them

Click [HERE<sup>4</sup>](#) to see the text message and Covid-19 Florey template  
from the patient's point of view

---

<sup>1</sup> <https://support accurx.com/en/articles/3879662-can-i-remotely-monitor-a-patient-s-covid-19-symptoms>

<sup>2</sup> <https://support accurx.com/en/articles/3548461-using-florey-to-gather-qof-data>

<sup>3</sup> <https://support accurx.com/en/articles/2678798-can-i-cancel-a-delayed-text-message-after-i-ve-scheduled-it>

<sup>4</sup> <https://youtu.be/RwVm7CpnaDU>